

blast off



eye	a	b	ed	dbl	temples
42	42.1	25.1	42.8	15	120mm
44	44.0	27.0	44.9	15	125mm

blue

brown

features

metal alloy
adjustable silicon unifit bridge
single action spring hinges
additional temple sizes available:
skull temples: 110, 115, 130mm
cable temples: 115, 120, 125, 130, 135cc

cotton candy



eye	a	b	ed	dbl	temples
36	35.9	24.9	36.1	15	115cc
38	37.9	26.8	37.9	15	120cc
40	39.8	28.8	39.9	15	125cc
42	42.0	30.9	42.1	15	130cc

blue

bronze

pink

features

metal alloy
cable temples standard
single action spring hinges
additional temple sizes available:
skull temples: 110, 115, 120, 125mm

ducky



eye	a	b	ed	dbl	temples
40	40.0	31.0	41.4	16	120cc
42	42.0	33.0	43.5	16	125cc
44	43.9	34.9	45.3	16	130cc

brown

lilac

features

acetate with silk screened ducky motif
bridge positions vary on different sizes for
developing nasal structure
cable temples standard
single action spring hinges
additional cable temples: 110, 115cc
additional skull temples: 110, 115, 120, 125mm

gold fish



eye	a	b	ed	dbl	temples
42	41.8	24.7	42.2	17	125mm
44	43.9	26.7	44.2	17	130mm

brown

gunmetal

features

metal alloy
gently bowed skull temples standard
single action spring hinges
additional temple sizes available:
skull temples: 115, 120, 135mm
cable temples: 115, 120, 125, 130, 135cc

lollipop



amber matte

gold matte

eye	a	b	ed	dbl	temples
34	33.7	25.9	33.9	17	115cc
36	35.7	28.3	35.8	17	120cc
38	37.7	30.1	37.9	17	125cc
40	40.1	32.3	40.1	17	130cc

features

metal alloy
cable temples standard
single action spring hinges
additional temple sizes available:
skull temples: 105, 110, 115, 120, 125mm

available July

luv bug



pink

violet

eye	a	b	ed	dbl	temples
38	38.0	23.4	39.0	15	120cc
40	39.8	24.7	41.1	15	125cc
42	42.0	25.9	43.2	15	130cc

features

metal alloy
cable temples standard
single action spring hinges
additional temple sizes available:
cable temples: 115cc
skull temples: 110, 115, 120, 125mm

marshmallow



brown

pink matte

eye	a	b	ed	dbl	temples
35	34.9	25.8	35.2	15	115cc
37	36.9	27.8	37.2	15	120cc
39	39.0	29.8	39.3	15	125cc
41	41.0	32.0	41.3	15	130cc

features

metal alloy
silicon unifix bridge
cable temples standard
single action spring hinges
additional temple sizes available:
skull temples: 110, 115, 120, 125mm

mary jane



pink


purple

eye	a	b	ed	dbl	temples
42	41.8	24.9	41.9	16	120mm
44	43.8	27.0	44.0	16	125mm

features

metal alloy
skull temples standard
single action spring hinges
additional temple sizes available:
skull temples: 110, 115, 130mm
cable temples: 115, 120, 125, 130, 135cc

peanut



bronze

pink

eye	a	b	ed	dbl	temples
35	34.8	23.8	35.2	16	115cc
37	36.9	24.9	37.5	16	120cc
39	38.8	25.8	39.1	16	125cc

features

metal alloy
 cable temples standard
 single action spring hinges
 additional temple sizes available:
 cable temples: 130cc
 skull temples: 105, 110, 115, 120, 125mm

teddy bear



blue

brown

pink

eye	a	b	ed	dbl	temples
36	35.2	31.4	36.3	15	115cc
39	38.2	34.4	39.3	15	120cc
42	41.2	37.4	42.3	15	120cc

features

acetate with silk screened teddy bear motif
 bridge positions vary on different sizes
 for developing nasal structure
 cable temples standard
 single action spring hinges
 additional temple sizes available:
 cable temples: 110, 125, 130cc
 skull temples: 100, 105, 110, 115, 120mm