

features and benefits **metal**



FEATURES AND BENEFITS IN CLEARVISION METAL STYLES

FEATURES

- Slotted (-) and/or (+) screw heads
- Special adjustable silicon unifit wing nose pad*
- 'D' shaped PVC click-on nose pads
- Hidden rim lock positions*
- Glued acetate temple sleeves
- Metal end piece for adjusting pantoscopic angle
- Correct tension of nylon cord holding grooved demo lenses
- Additional spooled "uncut" nylon cord with insertion ribbon for semi-rimless styles

HINGES

- Double action spring hinges*

DESIGN FEATURES


- Progressive friendly eye shapes
- Multiple temple lengths available
- Combination metal and acetate designs*
- Oversized ophthalmics*
- Age cross over

MATERIALS

- Metal alloy
- Combination metal and acetate designs

* featured in frame shown above: max gunmetal

features and benefits acetate



FEATURES AND BENEFITS IN CLEARVISION ACETATE STYLES

FEATURES

- Slotted (-) and/or (+) screw heads
- Flared nose pads for a perfect bridge fit*
- Metal end piece for adjusting pantoscopic angle

HINGES

- Double action spring hinges*

DESIGN FEATURES

- Progressive friendly eye shapes*
- Sculpted temples
- Special end piece designs*
- Combination metal and acetate designs*
- Wrapped end piece designs
- Accomodating base curve design
- Age cross over

MATERIALS

- Handmade acetate and combination designs
- Injected Cellulose Propionate plastics

*featured in frame shown above: jacqueline mauve