


features and benefits **metal**



FEATURES AND BENEFITS IN FISHER-PRICE METAL STYLES

FEATURES

- Custom molded 'Fisher-Price' logo PVC nose pads with rolled edges which are molded onto a screw-on pad arm attachment*
- Thicker eye rims and extra solder at bridge and temples for extra durability
- Large range of sizes for a custom fit
- Metal end piece for adjusting pantoscopic angle
- Strap bridge options in four sizes available for most metal frames*
- All frames undergo rigorous extensive testing at our Fisher-Price approved lab and surpass industry standards
- Single vision and multifocal Rx designed shapes

HINGES

- Heavy duty European spring hinges with single action for safety and durability*

DESIGN FEATURES

- Fisher-Price Eyewear was designed by pediatric specialists who understand the anatomy of children and their special fitting needs
- Multiple coiled comfort cable temple lengths available for all models*
- Multiple skull temple lengths available for all models*
- Skull temples feature exclusive "tip-up" safety system to prevent tips from falling off and a child swallowing
- Quality recognized branded frames for children

MATERIALS

- High nickel metal alloy materials that surpass government regulations

* featured in frame shown above: lollipop amber matte
inset image: mary jane purple

features and benefits acetate



HEAVY DUTY EUROPEAN SPRING HINGES WITH SINGLE ACTION FOR SAFETY AND DURABILITY

WIDE FLARED SADDLE BRIDGE WITH VARYING BRIDGE POSITION ON EACH SIZE APPROPRIATE FOR THE NASAL FIT FOR EACH DEVELOPING AGE

MULTIPLE CABLE AND SKULL TEMPLE LENGTHS AVAILABLE FOR ALL MODELS

STOCKY ACETATE MATERIAL FOR DURABILITY WITH SILK SCREENED TEDDY BEAR MOTIF



BOYS AND GIRLS COLORS AVAILABLE IN ALL MODELS

FEATURES AND BENEFITS IN FISHER-PRICE ACETATE STYLES

FEATURES

- Wide flared saddle bridge with varying bridge position on each size appropriate for the nasal fit for each developing age*
- Large range of sizes for a custom fit
- All frames undergo rigorous extensive testing at our Fisher-Price approved lab and surpass industry standards
- Single vision and multifocal Rx designed shapes

HINGES

- Heavy duty European spring hinges with single action for safety and durability*

DESIGN FEATURES

- Fisher-Price Eyewear was designed by pediatric specialists who understand the anatomy of children and their special fitting needs
- Multiple cable and skull temple lengths available for all models*
- Skull temples feature exclusive "tip-up" safety system to prevent tips from falling off and a child swallowing
- Boys and girls colors available*
- Silk screened motifs*
- Quality recognized branded frames for children

MATERIALS

- Stocky acetate material for durability*

* featured in frame shown above: teddy bear pink
insert image: teddy bear brown