

		Suggested Rx Range PD 60 - 64mm to -3.50 cylinder
Model	Style	
1109	ARTHUR	+5.50 to -5.00
1153	CANYON	8 BASE: +3.50 to -3.50
1092	CAPER	8 BASE: +3.50 to -4.00
1154	CLIVE	+5.00 to -4.00
1106	CONRAD	+5.00 to -5.00
1164	COVE	+5.50 to -5.50
1027	CRAWLER	+3.00 to -4.00
1101	DAPHNE	+5.00 to -6.00 ^
1084	DECOY	+3.50 to -3.50
1169	DESCEND/A	8 BASE: +2.50 to -3.00 6 BASE: -3.25 to -4.00 †
4060	DESCEND/E	8 BASE: +2.50 to -3.00 6 BASE: -3.25 to -4.00 †
4059	DESCEND/N	8 BASE: +2.50 to -3.00 6 BASE: -3.25 to -4.00 †
1168	DESCEND/O	8 BASE: +2.50 to -3.00 6 BASE: -3.25 to -4.00 †
1070	DESCEND XL	8 BASE: +2.50 to -3.00 6 BASE: -3.25 to -4.00 †
1070	DESCEND/Z	AVAILABLE AS PLANO ONLY
1165	DESERT	+5.00 to -5.00
1127	DEXTER	+5.00 to -4.00
1097	ESPEN	8 BASE: +3.00 to -3.00
1112	FINLEY	+5.00 to -6.00
1096	FORGE	+5.00 to -5.00
1145	HANK	+3.50 to -3.50
1082	HARBOR	+6.00 to -6.00
1108	HARRISON	+5.50 to -5.00
4071	HARNESS	+2.00 to -3.50
1019	HOLSBY	+2.00 to -3.50
1111	JASPER	8 BASE: +5.00 to -4.00
1069	LEIGH	+4.00 to -5.00
1098	MAVERICK	8 BASE: +3.00 to -3.00
1080	MAXIE	+6.00 to -6.00
1163	METRO	+4.50 to -4.50
1139	PROSPER	8 BASE: +2.50 to -3.00
1135	QUINN	+5.50 to -5.00
1011	RACONTEUR	+2.50 to -4.00
1137	REBEL	8 BASE: +2.50 to -3.00
1014	RELAY	+2.50 to -4.00
1156	RELAY PETITE	+5.00 to -5.00
1143	RILEY	+5.00 to -6.00
1136	ROSE	+5.00 to -6.00 ^^
1102	SAMMY	+5.00 to -6.00 ^
1161	SIERRA	+5.50 to -6.00
1038	SKYLAR	+2.00 to -3.00
1050	SLATER	+3.00 to -4.50
1081	SPARK	+6.00 to -6.00
1116	SUMMIT	8 BASE: +3.50 to -3.50
1138	SURGE	8 BASE: +2.50 to -3.00
1083	TARQUIN	8 BASE: +3.50 to -3.50
1079	TATE	8 BASE: +3.50 to -3.50
1104	TAYLOR	+5.00 to -6.00
1129	WATSON	+5.50 to -5.00
1054	ZINGER	-4.50 to +3.00

* Rx range listed is for sphere and cylinder combined (minus cylinder form.) ** Due to multiple variables, certain Rx/PD combinations may not be possible

*** Ranges listed are for polarized polycarbonate. NXT may vary.

† Descend series note: Up to -3.00 will be 8 base. From -3.25 to -4.00 will be on 6 base with slightly adjusted drill points.

^ Minimum PD 55, ^^ Minimum PD 58